

IHE Actors and Transactions

Profile	Actor	Transaction	PB Supported
Scheduled Workflow (SWF)	Image Manager	Procedure Schedule	v3.2
		Modality PPS in Progress	v6.0
		Modality PPS Completed	v6.0
		Storage Commitment	v5.0
		Image Availability Query	v3.1
		Patient Update	v4.0
		Creator PPS in Progress	v6.0
	Image Archive	Creator PPS Completed	v6.0
		Modality Image Stored	v3.0
		Query Images	v3.1
		Retrieve Images	v3.1
		Creator Image Stored	v3.0
		Evidence Documents Stored	v3.0
Patient Information Reconciliation (PIR)	Image Manager	Procedure Schedule	v3.2
		Modality PPS in Progress	v6.0
		Modality PPS Completed	v6.0
		Storage Commitment	v5.0
		Image Availability Query	v3.1
		Patient Update	v4.0
		Procedure Update	v4.0
	Image Archive	Modality Image Stored	v3.0
Consistent Presetation of Images (CPI)	Image Manager	Modality PPS in Progress	v6.0
		Modality PPS Completed	v6.0
		Storage Commitment	v5.0
		Creator PPS in Progress	v6.0
		Creator PPS Completed	v6.0
	Image Archive	Modality Image Stored	v3.0
		Modality Pres State Stored	v4.1
		Query Images	v3.1
		Query Pres State	v4.1
		Retrieve Images	v3.1
		Retrieve Pres State	v4.1
		Creator Image Stored	v3.0
		Creator Pres State Stored	v4.1
	Print Composer	Print Request with LUT	v3.0
Presentation of Grouped Procedures (PGP)	Image Manager	Procedure Scheduled	v3.2
		Modality PPS in Progress	v6.0
		Modality PPS Completed	v6.0
		Storage Commitment	v5.0
		Creator PPS in Progress	v6.0
		Creator PPS Completed	v6.0
		Image Availability Query	v3.1
	Image Archive	Modality Image Stored	v3.0
		Modality Pres State Stored	v4.1
Access to Radiology Information (ARI)	Image Archive	Query Images	v3.1
		Query Pres State	v4.1
		Retrieve Images	v3.1
		Retrieve Pres State	v4.1

IHE Actors and Transactions

	Report Repository	Query Reports	v4.1
		Retrieve Reports	v4.1
Key Image Note (KIN)	Image Manager	Modality PPS in Progress	v6.0
		Modality PPS Completed	v6.0
		Storage Commitment	v5.0
		Creator PPS in Progress	v6.0
		Creator PPS Completed	v6.0
	Image Archive	Query Images	v3.1
		Query Pres State	v4.1
		Retrieve Images	v3.1
		Retrieve Pres State	v4.1
		Store Key Image Note	
		Query Key Image Note	
		Retrieve Key Image Note	
Simple Image and Numeric Report (SINR)	Report Creator	Report Submission	v4.0
	Report Manager	Report Submission	v4.0
		Report Issuing	v4.1
		Structured Report Export	v4.0
	Report Repository	Report Issuing	v4.1
		Query Reports	v4.1
		Retrieve Reports	v4.1
Basic Security (SEC)	Image Manager	Record Order Event	
	Image Archive	Record Order Event	
	Print Composer	Record Order Event	
	Report Creator	Record Order Event	
	Report Manager	Record Order Event	
	Report Repository	Record Order Event	
Charge Posting Transaction Diagram (CHG)	N/A		
Post-Processing Workflow (PWF)	Image Manager/ Image Archive	Storage Commitment	v5.0
		Images Availability Query	v3.1
		Patient Update	v4.0
		Query Images	v3.1
		Retrieve Images	v3.1
		Creator Images Stored	v3.0
		Evidence Documents Stored	v3.0
		Query Evidence Documents	v3.1
		Retrieve Evidence Documents	v3.1
		Perf Work Status Upd (Send)	
Reporting Workflow (RWF)	TBD		
Evidence Document (ED)	TBD		
Portable Data for Imaging (PDI)	TBD		
NM Image (NM)	TBD		